

model # 481015

LS&S, LLC  
145 River Rock Dr  
Buffalo, NY 14207-2172

**TAYLOR**

PRECISION  
**TECH**<sup>TM</sup>

Lithium  
Electronic Scale  
with Voice Feature

## Instruction Manual

Congratulations on your purchase of a Taylor® Lithium Electronic Scale with Voice Feature. Your Scale is an example of superior design and craftsmanship. Please read this instruction manual carefully before use. Keep these instructions handy for future reference.

### General Operation

---

1. There is a red strip under the lithium battery cover to prevent battery drain in shipping. The strip can be pulled out without removing the battery cover. Also, the 9 volt battery under the cover that has "For Voice" printed on it has to be removed to connect it. Open the cover by pushing in the direction of the arrow, lift and remove the cover. The battery has an adhesive strip over the terminals that must be removed before snapping the connector on the battery. Follow the polarity symbols marked on the outside of the cover. Place the battery back in the compartment and be careful not to pinch the wires as you snap the cover in place. Please do this before continuing with these instructions.
2. Select pound, kilogram or stone-pound by the slide switch on the bottom of the scale.
3. Select language of English, French, or Spanish by pressing the button on the bottom of the scale. Sound will output the corresponding language once the language is selected. You can also select a volume level of high, low or off by using the slide switch on the bottom of the scale.

**"Note: English, French and Spanish can be selected for kilogram units. English and Spanish can be selected for pounds units. Stone units can be heard in English only."**

## General Operation (continued)

---

4. *The scale needs to be initialized at the first time of use or after any battery replacement: Press lightly on the scale platform. The display will show zero. Then the scale automatically turns off. Your scale is now ready for use.*
5. *With the automatic on function, simply step onto the scale, and you will hear "Hello". After a few seconds it will show your weight and you will hear "Your weight is -----".*

## Trouble Shooting

---

1. The message "Err0" is displayed with sound "Error, please try again." If this happens simply step off the scale, and repeat the weight measurement procedure.
2. In the event that you are moving while you are standing on the platform, the scale will display the message "Err1" with sound "Error, please try again!" If this happens, step off the scale, then repeat the measurement procedure ensuring that you keep perfectly still.
3. **OVERLOAD WARNING:** Remove the weight immediately when the scale displays "Err2" with sound "Overload", Otherwise, permanent damage to the scale will occur.

## Battery Replacement

---

1. *Replace lithium batteries when **BATT** is displayed or replace 9 volt battery when the voice quality distorts or becomes too soft to hear.*
2. *Use two CR2032 Lithium batteries for weighing and one 9 volt alkaline battery for voice.*
3. *Insert batteries into the appropriate compartments on the bottom of the scale as indicated by the polarity symbols(+ and -) as marked on the outside of the battery compartments. After replacing any battery follow step 4 under "General Operations" to initialize the weight reading.*
4. *Remove batteries if the scale will not be used for a prolonged period of time.*

# Reaching Target Weight

Maintaining proper weight is essential to healthy living. Obesity has been linked to a higher risk of heart disease, diabetes and some forms of cancer. The Body Mass Index (BMI) is the most widely accepted measurement of weight according to health professionals. Use the chart below to determine your own BMI by matching your height in the left hand column with your weight in the center. Anyone, male or female, with a BMI of 25 to 29 is considered over weight and a body mass index of 30 is considered obese. Remember to eat right and exercise.

## Body Mass Index

Height (Inches)	Body Weight (pounds)						Overweight						Obese					
	91	96	100	105	110	115	119	124	129	134	138	143	148	153	158	162	167	
58	91	96	100	105	110	115	119	124	129	134	138	143	148	153	158	162	167	
59	94	99	104	109	114	119	124	128	133	138	143	148	153	158	163	168	173	
60	97	102	107	112	118	123	128	133	138	143	148	153	158	163	168	174	179	
61	100	106	111	116	122	127	132	137	143	148	153	158	164	169	174	180	185	
62	104	109	115	120	126	131	136	142	147	153	158	164	169	175	180	186	191	
63	107	113	118	124	130	135	141	146	152	158	163	169	175	180	186	191	197	
64	110	116	122	128	134	140	145	151	157	163	169	174	180	186	192	197	204	
65	114	120	126	132	138	144	150	156	162	168	174	180	186	192	198	204	210	
66	118	124	130	136	142	148	155	161	167	173	179	186	192	198	204	210	216	
67	121	127	134	140	146	153	159	166	172	178	185	191	198	204	211	217	223	
68	125	131	138	144	151	158	164	171	177	184	190	197	203	210	216	223	230	
69	128	135	142	149	155	162	169	176	182	189	196	203	209	216	223	230	236	
70	132	139	146	153	160	167	174	181	188	195	202	209	216	222	229	236	243	
71	136	143	150	157	165	172	179	186	193	200	208	215	222	229	236	243	250	
72	140	147	154	162	169	177	184	191	199	206	213	221	228	235	242	250	258	
73	144	151	159	166	174	182	189	197	204	212	219	227	235	242	250	257	265	
74	148	155	163	171	179	186	194	202	210	218	225	233	241	249	256	264	272	
75	152	160	168	176	184	192	200	208	216	224	232	240	248	256	264	272	279	
76	156	164	172	180	189	197	205	213	221	230	238	246	254	263	271	279	287	
	<b>19</b>	<b>20</b>	<b>21</b>	<b>22</b>	<b>23</b>	<b>24</b>	<b>25</b>	<b>26</b>	<b>27</b>	<b>28</b>	<b>29</b>	<b>30</b>	<b>31</b>	<b>32</b>	<b>33</b>	<b>34</b>	<b>35</b>	

## ***Precautions***

---

1. *Clean exposed parts with a soft, slightly damp cloth. To remove stains, use a mild soap. Never use detergents, excess water, treated cloths, harsh cleaning agents, or sprays. Do not attempt to lubricate or open scale casing as this will void your warranty.*
2. *Although your scale is designed to be maintenance-free, it contains sensitive electronic parts. Avoid rough treatment. Do not drop, kick, or jump on it. Treat it with care to ensure the best performance.*
3. *This scale is designed to read weights between 25 and 330 lbs./12 and 150 kgs and on special models 1 and 23 stone. To prevent damage do not place more than 330 lbs./150 kgs./23 stone on it.*
4. *Do not store scale where you store cleaning chemicals. The vapors from some household products may affect the electronic components of your scale. Do not store scale on its side.*
5. *This scale is an extremely sensitive weighing device. In order to prevent run down of the battery, do not place anything on the scale.*

## ***Lifetime Warranty***

---

This scale is warranted against defects in materials of workmanship for the life of the original purchaser from date of retail purchase. It does not cover damages or wear resulting from accident, misuse, abuse, commercial use, or unauthorized adjustment and/or repair.

Do not return to retailer. Should this scale require service (or replacement at our option) while under warranty, please pack the item in the original packaging and return it prepaid, along with store receipt showing date of purchase and a note explaining reason for return to:

**TAYLOR.**    **PRECISION PRODUCTS**

---

**2220 Entrada Del Sol  
Las Cruces, New Mexico 88001  
www.taylorusa.com**

There are no express warranties except as listed above. This warranty gives you specific legal rights, and you may have other rights which vary from state to state.

©2003 Taylor Precision Group, LLC and its affiliated companies, all right reserved. Taylor® is registered trademarks of Taylor Precision Group, LLC and its affiliated companies. Precision Tech is a trademark of Taylor Precision Group, LLC and its affiliated companies. All rights reserved.

**Made to our exact specifications in China.**